


# Record Calcium & Contractility from Intact Cardiomyocytes



Complete real-time,  
turnkey system

Fast, quantitative calcium,  
sarcomere shortening,  
and cell length


Intuitive software for  
data acquisition and  
transient analysis

## IonOptix Calcium & Contractility System

The IonOptix Calcium & Contractility System offers simultaneous acquisition of fluorescence photometry with digital cell geometry measurements.

All components are designed and engineered to work as part of a complete and integrated workstation for precise, synchronous data acquisition.

Systems are fully installed on-site with extensive training so researchers can be ready to collect data within days of delivery. And customers can rest easy with unlimited support for the system lifetime.





Typical C&C System with MyoPacer cell stimulator, HyperSwitch, and IonWizard acquisition module

## Acquisition Features

- **Maximum sampling rates** of 1000 Hz for cell length and sarcomere data; 1000 Hz for single excitation calcium; 250 Hz for ratiometric dual excitation calcium
- Separate acquisition time periods defined as '**epochs**' that can each include independent sampling rates and analog data. Seamlessly switch between epochs automatically or with a trigger
- **Cell length** algorithm uses image intensity or derivative of image intensity to detect left and right edges
- **Sarcomere spacing** algorithm uses fast Fourier transform for reliable, reproducible real-time measurements
- Analog monitoring voltages and **synchronization signals** (e.g. pacing sync) collected simultaneously

## Analysis Parameters

- **Departure/return velocity (d/dt):** maximum/minimum velocity reached on the rising/falling phases of the transient
- **Peak height:** can be expressed as absolute peak height or % of the baseline (e.g. % shortening for length data)
- **Time to % peak or time to % baseline:** time required for transient to reach a specified level on the rising and falling phases (e.g. TP90 and TR90)
- **Exponential fit:** single exponential curve fit and associated time constant ( $\tau$ ), useful for characterizing the rate of calcium re-uptake
- **Area under the curve:** useful for measuring total calcium exchanged
- **Event averaging:** average data to remove noise from repetitive signals



## Select Publications

Cardiac resynchronization therapy restores sympathovagal balance in the failing heart by differential remodeling of cholinergic signaling.

DeMazumder, Deeptankar, et al. Circulation research 116.10 (2015): 1691-1699. DOI: 10.1161/CIRCRESAHA.116.305268

Mechanochemotransduction during cardiomyocyte contraction is mediated by localized nitric oxide signaling.

Jian, Zhong, et al. Science signaling 7.317 (2014): ra27. DOI: 10.1126/scisignal.2005046