



MURINE VENTRICULAR CARDIOMYOCYTE ISOLATION, CALCIUM ADAPTATION, AND FURA-2/AM DYE LOADING PROTOCOL

Isolation Solutions:

1. Prepare stock solutions: 0-Ca²⁺ Physiologic Saline Solution (PSS), MEM Base, and 2mM-Ca²⁺ PSS (see solution supplement)
2. Prepare enzyme aliquot (see solution supplement)

Immediately Prior to Isolation:

3. Make isolation solutions (see solution supplement)
 - a. Wait to add BSA to solution 3 until immediately prior to anesthetizing the rodent – pre-measure and store at 5° C.
4. Check oxygen tanks
5. Rinse DI water through isolation rig
 - a. Cover top reservoir until use
6. Check vacuum flasks
7. Check water level in recirculation pump
 - a. Check bath temp
8. Boil 1L Millipore water to clean rig after isolation
9. Turn on water bath, recirculation pump, and temperature controllers
10. Add solution 1 to top reservoir and turn on oxygen to outgas
11. Perfuse solution 1 through tubing, removing all bubbles
12. Tie double knot sutures onto cannula

Cannulation and Isolation Protocol:

13. After sufficiently anesthetizing rodent, excise the heart per lab protocol (varies by ACUC)
14. Place heart in cold 0-Ca²⁺ PSS, removing excess tissue as needed
15. Cannulate aorta and secure with pre-tied suture; begin solution 1 perfusion
16. Perfuse solution 1 ~2 minutes, ensuring heart fully submerged in solution
17. During this time, add enzyme to solution 2
18. Perfuse solution 2 and set enzymatic digestion timer to coincide with solution arrival at heart
19. Once solution 2 is in the heart, drain the bath of all excess solution 1 via transfer pipette to ensure heart is bathed only in solution 2
20. Solution 2 perfusion time is ~6-12 minutes (highly condition- and tissue-dependent). Barely visible cellular detritus should be apparent after sufficient digestion. Tissue should be orange in color and soft but not slimy.
21. Fill weigh boats with solution 3 (stop solution with BSA)
22. After sufficient digestion, carefully remove heart from cannula and submerge in solution 3
 - a. Remove atria, excess fat, and right ventricle (with or without septum) via fine surgical scissors
 - b. Move trimmed left ventricle to separate weigh boat and make 5-6 fine cuts to help dissociate cells



- c. Gently agitate with fine forceps for ~5-10 minutes to further dissociate cells
 - d. Check cells under microscope for live/dead ratio
 - e. Remove heart by scraping gently up the side of the weigh boat to further remove cells
23. Cells are filtered into 15ml falcon tube through 100-200 μm mesh filter
- a. A few drops of solution 3 to pre-wet will break static seal on mesh
24. Use an angled surface or box to keep tube at $\sim 45^\circ$ angle during all subsequent procedures
- a. Failure to allow pelleting at an angle can result in cell clumping and sarcolemmal damage (followed by calcium overload)
25. Allow cells to pellet ~ 10 minutes, remove supernatant, and replace with fresh BSA/Stop solution to bring cells to desired density (solution 3)
26. Excess tissue can be discarded appropriately

Following Isolation (this can be done while plating/dye loading cells):

27. Once heart is removed from rig: oxygen, water bath, and recirculation pump can be turned off.
28. Wash rig (make sure vacuum dish is open):
- a. Flush with $\sim 500\text{ml}$ of boiling Millipore water through rig
 - b. Flush with $\sim 500\text{ml}$ of 70% ethanol
 - c. Flush again with $\sim 500\text{ml}$ of boiling Millipore water
 - d. Flush a final time with Millipore water to ensure all ethanol is removed
29. After all solution has been flushed, turn off hotplate and vacuum.

Assessing Cardiomyocyte Isolation Quality:

30. Assess newly isolated cardiomyocytes for live:dead ratio
- a. $>60\%$ represents experimentally viable cells
 - b. This percentage is model- and protocol-dependent as isolations in disease models (e.g. aging, muscular dystrophy, western diet, etc.) may have significantly lower yields
31. Visually inspect isolated cardiomyocytes under bright field for:
- a. Distinct sarcomeric striations
 - b. Lack of sarcolemmal “blebbing” (bubble-like visual indicators of compromised sarcolemma integrity)
 - c. Lack of non-stimulated cardiomyocyte contraction (spontaneous contraction), indicative of compromised sarcolemmal integrity
 - d. Crisply defined cellular borders (not “chewed up” or clumped together with other cells, indicative of over- and under-digestion, respectively)
 - e. If time permits, resting sarcomere lengths $>1.8 \mu\text{m}$ indicate experimentally viable cells
32. Critically, ensure visual markers of cellular integrity are intact following adaptation to physiologic calcium concentration and cells are not spontaneously contracting.
- a. Physiologic calcium adapted cardiomyocyte functional integrity is the *critical* hallmark of a good isolation

Adaptation and Plating

33. Re-suspend freshly isolated cardiomyocytes in fresh BSA/STOP solution after removing old BSA
34. Step-wise over 20 minutes, adapt cells from nominal to $250 \mu\text{M Ca}^{2+}$
- a. Suggested interval: 5 minutes



- b. Lightly re-suspend cells prior to additions using large-bore transfer pipette
35. Allow cells to pellet, remove supernatant, and re-suspend in 500 μM Ca^{2+} PSS
- a. *NOTE:* Cells can be further adapted and loaded either plated on laminin-coated coverslips or suspended in falcon tubes. Read further for plating instructions:
 - i. For coverslips, coat with $\sim 3\text{-}5$ μL 1M laminin and smear evenly in appropriate area.
 - ii. Once completely dry, transfer pipette sufficient suspended cardiomyocytes to cover laminin-coated area. To minimize dessication, place moistened pieces of kimwipe in petri dishes containing slides and cover with lid.
 - iii. Allow 25-30 minutes for cardiomyocyte attachment.
 - iv. Below, ignore pelleting instructions and simply pipette off excess solution between solution changes, replacing with indicated solution
36. After 5 minutes, adapt cells to 750 μM Ca^{2+}
37. Allow cells to pellet, remove supernatant, and re-suspend in 1 mM Ca^{2+} PSS

Fura-2/AM Dye Loading

38. Dye load isolated cardiomyocytes in 2-5 μM fura-2/AM in DMSO (empirically derived)
39. Allow cells to load for 20-25 minutes (empirically derived)
40. Wash (de-esterification of acetoxymethyl ester) with 2 mM Ca^{2+} for ~ 30 minutes
- a. Ensure near-complete removal of fura-2/AM solution prior to washing
41. Cells are ready for experimentation

Adam B. Veteto, Ph.D
Applications Scientist
IonOptix LLC
Office: (617) 696-7335
Mobile: (781) 366-8815
Email: adam@ionoptix.com
Web: www.ionoptix.com

Michelle D. Lambert, BS
Research Specialist
University of Missouri – Columbia

Timothy L. Domeier, Ph.D
Associate Professor
University of Missouri – Columbia