



# MCS300

# MyoCam-S3

User Manual

*This page is intentionally left blank.*

## **Copyright 2022 IonOptix LLC**

All rights reserved. No parts of this work may be reproduced in any form or by any means - graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems - without the written permission of the publisher.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.

Document date: November 8, 2022

Printed: April 17, 2026 in Westwood, Massachusetts, USA

### **IonOptix LLC**

396 University Ave

Westwood, MA 02090

USA

<https://www.ionoptix.com/>

[info@ionoptix.com](mailto:info@ionoptix.com)

US: +1 617 696 7335

EU: + +31 317 728 104

### **Research Use Only**

This product is intended for research purposes only. It is not certified for clinical applications (including diagnostic purposes). Use of this product in uncertified applications is in violation of FDA regulations.

# Table of Contents

Safety Precautions	1
Statement of Conformity	1
Intended Use	1
Disposal	2
Warranty	2
Technical Specifications	2
Overview	3
Features	3
Connections	3
USB3	3
IN (Timing Trigger Input)	3
Out (Frame Active Signal)	3
Installation	4
How to Use	4
Cleaning and Maintenance	4
Cleaning	4
Maintenance	4
Storage and Transport	4
Storage	4
Transport	4
Faults and Troubleshooting	4
Troubleshooting	4
Repairing	5
Accessories	5
Index	6

## Safety Precautions



To ensure reliability and safety, we recommend using the instrument according to the guidelines described in this manual. All safety and operating instructions should be read and understood before use.

Electric shock risk! Do not remove instrument cover. There are no user-serviceable parts inside. Refer servicing to IonOptix technical support ([support@ionoptix.com](mailto:support@ionoptix.com)). Static electricity can damage electronic components. Take care to discharge yourself before handling the device.

Please adhere to the following safety guidelines:

- Unplug the instrument from the wall before cleaning.
- Keep the instrument free from moisture, water, and dust.
- Do not place instrument on unstable surface and do not drop the device.
- If ventilation slots are present on instrument panels, keep clear of obstructions.
- All cables and power supply cords should be routed so as to not cause a hazard. Use only supplied or approved/recommended cables.
- Ensure instrument is not situated near a heat source.

## Statement of Conformity



This instrument carries the CE mark and meets the appropriate EU directives.

## Intended Use



This product is intended for **research purposes only**. It is not certified for clinical applications (including diagnostic purposes). Use of this product in uncertified applications is in violation of FDA regulations.

IonOptix products are intended to be installed, used, and operated under the supervision of an appropriately qualified life-science researcher.

## Disposal



Forward to recycling center or return to manufacturer. Unwanted equipment bearing the Waste Electrical and Electronic Equipment (WEEE) Directive symbol requires separate waste collection. For a product labeled with this symbol, either forward to a recycling center or contact your nearest IonOptix representative for methods of disposal at the end of its working life.

## Warranty

This IonOptix instrument is warranted for a period of one (1) year from delivery under normal use and conditions, and will operate in material conformance with applicable IonOptix specifications or standards.

This warranty does not apply to any defect caused by failure to provide a suitable operating environment, any other excluded use, or any other abuse, misuse, or neglect of the instrument.

## Technical Specifications

Item	Specification
Dimensions	82 (L) x 45 (W) x 45 (H) mm
Weight	175 g
Resolution	2.4 MPixel (1936 x 1216)
Connection	USB3
IN / OUT connector physical standard	BNC
IN / OUT connector electrical standard	5V TTL/CMOS
Power	USB bus powered (5V)
Sensor	1/1.2" CMOS (5.86 $\mu$ m pixels)
Operating temperature	10 - 40 °C
Air humidity	10 – 90 % RH, non-condensing
Pollution degree	2
Exposure Time	10 $\mu$ s – 20 s
Frame Rate	161.4 Hz (subfield frame rates > 1 kHz)
ADC Resolution	12 bit
SNR (max)	45.1 dB (EMVA1288)
Dynamic range	66.4 dB (EMVA1288)
Programmable gain and offset	USB3 Vision industry standard
Lens Mount	C-mount

## Overview

This manual describes how to setup and use the MyoCam-S3 camera in your application.

The IonOptix MyoCam-S3™ is our third generation MyoCam. Featuring a faster, higher resolution CMOS sensor, it is uniquely designed for use with IonOptix Contractility Systems where high frame rates and sub-micron resolution are required to capture myocyte shortening.

The MyoCam-S3™ is an all-digital, variable field rate camera that utilizes the USB3 Vision standard for higher data transfer rates over its USB2 predecessor while retaining an easy-to-use plug and play interface.


Power and data occur over the same lockable USB3 cable, enabling simple, secure setup and connectivity with better performance and a lower performance overhead.

The MyoCam-S3™ gives you complete control of all aspects of video acquisition to deliver the optimum combination of temporal and spatial resolution, making it the ideal choice for superior myocyte data acquisition.

## Features

- Single, lockable USB3 cable for power and data transfer
- Compact form factor
- Programmable integration time to stop fast movement or increase camera sensitivity
- Programmable frame period and external trigger to synchronize with other recording devices

## Connections



### USB3

USB3 connection to computer. Use the recommended cable with jack screws to prevent inadvertent disconnections.

### IN (Timing Trigger Input)

This is a 5 Volt TTL/CMOS logic level input used to synchronize the camera timestamps with the rest of the fluorescence system interface data streams.

### Out (Frame Active Signal)

This is a 5 Volt TTL/CMOS logic level output that is high when the sensor is actively integrating photons. This can be used as a trigger or for debugging purposes

## Installation



**Attention:** Always use the supplied or recommended cables and mounting accessories.

- Install camera with the C-Mount into microscope or CFA.
- Connect USB cable (USB3) to computer.
- Connect IN to Start Out of IonOptix Fluorescence System Interface (FSI). FSI must be FSI800 or FSI700 with FPGA circuit (manufactured in 2014 or later).

## How to Use

There are no user controls on the camera itself. All options are controlled by software. See the IonWizard Acquisition Manual for software setup and control options.

## Cleaning and Maintenance

### Cleaning

- Use a soft cloth and a mild cleaner (e.g. window cleaner) to clean the outside of the camera as necessary.
- Do not use solvents like acetone or thinner.
- Do not allow device to become saturated with cleaning solution.
- Do not attempt to clean the image sensor!

### Maintenance

- No maintenance other than cleaning is necessary.
- In case of malfunctions contact the manufacturer

## Storage and Transport

### Storage

- When not in use, keep the camera in the box it was delivered in.
- Store in a dry and clean environment

### Transport

- For transport, the original box is preferred.

## Faults and Troubleshooting

### Troubleshooting

- If the image stream does not start, check USB cable connection and FSI Start-out to Camera In connection.

## Repairing

- If the MyoCam-S3™ is damaged, stop using it and disconnect all cabling.
- The MyoCam-S3™ cannot be repaired by the user, and there are no user-serviceable parts inside.
- For repair or inspection, contact the manufacturer.

## Accessories

Part Number	Description
Z-C-USB-3.0-MAB	USB 3.0 A/M to Micro B/M with Dual Screw Lock Cable, 3m (10 ft.)
Z-C-BNC-MM-6FT	6 ft, male to male bnc cable
Z-MC-CF-MHF-S3	Myohandle fork for MCS300
Z-MC-CM/UC	C-mount/UC adapter
Z-OS-1mm/100Div	Stage Micrometer

# Index

## A

Accessories, 5

## C

Cleaning tips, 4  
Configuration, 4  
Connections, 3

## D

Declarations, 1  
Disposal, 2

## I

Intended Use, 1  
Introduction, 3

## S

Safety, 1  
Specifications, 2  
Storage, 4  
Synchronization, 3

## T

Troubleshooting, 4

## U

Usage, 4

## W

Warranty, 2